

Specialists in endoscope reprocessing
for more than 30 years!

EndoHigh® GTA
High level disinfectant for the
WASSENBURG® Endoscope
washer-disinfectors

EndoHigh® Detergent
The optimum detergent for the
WASSENBURG® Endoscope
washer-disinfectors

+ **Fast process: 28 minutes!***

New products with excellent benefits:
fast and effective.

+ **Safety first**

Fully tested detergent and disinfectant
with excellent efficacy, according to
ISO 15883.

+ **Brilliant combination**

Specially designed for use in combination
with WASSENBURG® Endoscope washer-
disinfectors.



WASSENBURG
medical

PROVEN PROCESS

EndoHigh® GTA has been specially developed for high level disinfection of endoscopes reprocessed in the WASSENBURG® WD440 and WASSENBURG® WD440 PT (pass through) Endoscope washer-disinfectors. EndoHigh® GTA complies with the requirements for disinfecting flexible endoscopes according to ISO 15883-4. The automated reprocessing process with EndoHigh® Detergent and EndoHigh® GTA meets the requirements of ISO 15883-4 regarding > 9 log microbiological reduction for the complete process.

ENDOHIGH® GTA SPECIFICATIONS

Active ingredient	Glutaral 10.5%
Dosing concentration	1.0%
Efficacy claim	5 min at 55-60°C
Dilution required	automatically performed in the WASSENBURG® washer-disinfector. The EndoHigh® GTA is for single use and no activator is required.

ENDOHIGH® DETERGENT SPECIFICATIONS

Active ingredient	Optimized blend of detergents and solubilizers
Dosing concentration	0.5%
Efficacy claim	4 min at 35-40°C
Dilution required	automatically performed in the WASSENBURG® washer-disinfector

Using EndoHigh® Detergent for manual pre-cleaning is preferred.



Only 28 minutes!*

Safe and reliable reprocessing only takes 28 minutes per endoscope.



Extra safe for healthcare professionals

The low GTA level (10.5%) and the mild odour mean, safe handling for the operator.



Excellent process, strong evidence

Fully tested with the latest test methods described in the ISO 15883 standard.



Great endoscope compatibility

The extended endoscope material compatibility testing give a wide endoscope compatibility.

* Exact process times depend on the machine model and program selected.

Please contact your Wassenburg supplier for more information on the EndoHigh® GTA and EndoHigh® Detergent.

© Copyright 2015. Subject to alterations.