

Specialists in endoscope reprocessing
for more than 30 years!

WASSENBURG® DRY300

Drying and storage cabinet



- + **Care for your endoscope**
Safe & controlled drying and conditioned storage for up to 10 flexible endoscopes.
- + **Always endoscopes ready for use**
Prevents recontamination and bacterial growth after washing & disinfection.
- + **Suitable for every logistic situation**
Available as pass through and non-pass through systems.
- + **Know where your endoscopes are**
One traceability system for both washing & disinfection and drying & conditioning.
- + **Pass through**
The DRY300 series is also available as a pass through system.

Also available with
integral compressor

OPTIMUM CONTROL

- + Endoscope storage time control
- + Temperature control
- + Filter condition control
- + Air circulation control
- + Independent drying positions
- + Bacteria free air through channels

DIMENSIONS AND WEIGHT


Width x Depth x Height: 1300 x 450 x 2070 mm
 Width with doors open (180°): 2000 mm
 Depth with doors open (90°): 960 mm
 Weight when empty (approx.): 280 kg
 Clearance required for access during service above the machine (approx.): 50 cm
 Right and left sides of the machine: the systems can be placed next to each other

AIR SUPPLY

Compressed air pressure: min. 2.0 bar, max. 8.0 bar
 Compressed air quality: oil-free, filtered, non-condensing

ELECTRICAL CONNECTION

Voltage: 220 – 240 V (frequency 50Hz)
 Type of electrical connection: plug with ground connection



Please contact us for the full installation specifications of the WASSENBURG® DRY300 and DRY300 D drying cabinets.
 The WASSENBURG® DRY300 and DRY300 D are CE-marked. © Copyright 2015. Subject to alterations.

Wassenburg Medical B.V.
 Edisonring 9
 6669 NA DODEWAARD
 THE NETHERLANDS

T: +31 (0) 488 700 500
 F: +31 (0) 488 453 685
 info@wassenburgmedical.nl
 www.wassenburgmedical.com