

Specialists in endoscope reprocessing since 1984!



## ENDOFENDER®

Offering a safe process in  
protecting expensive equipment.

## A 3 IN 1 SOLUTION

The EndoFender® developed by Wassenburg Medical features an innovative and safe technology that protects the delicate components of the endoscope and keeps the valves as a unique set with the parent endoscope during the entire reprocessing cycle. Handling this single use light weight device is very easy; the ergonomic design makes it possible to simply open and close it with one hand.



## HOW DOES IT WORK?

The three steps of the EndoFender® system cover the entire reprocessing cycle. The green sliding button indicates at which step the endoscope is within the reprocessing cycle.



- I The EndoFender® enters the reprocessing cycle right after the valves being manually brushed according to endoscope manufacturer's instructions. Valves are stored in opened position to be cleaned and disinfected inside EWD.
- II Containing endoscope and valves after being cleaned and disinfected inside EWD.
- III Used valves stay inside endoscope after procedure. Endoscope returns to reprocessing room.

## ENDOFENDER® PROTECTION DURING TRANSPORT AND STORAGE

Protect your expensive equipment with this tip protector; the EndoFender® not only protects the distal end but also the entire Bending Section Rubber, saving you time and money on repairs.

Valves can be placed in an open state in any of the two positions inside the EndoFender®. Distal caps can be stored inside the EndoFender® as well.



### 3 IN 1 SOLUTION

- Protecting the tip and Bending Section Rubber during drying, storage and transport.
- Keeping valves as a unique set with parent endoscope; better traceability.
- Valves are opened during cleaning & disinfection inside washer-disinfector.

## COMPLIANCY

EndoFender® keeps you compliant with guidelines stating that valves should be kept as a unique set with their corresponding endoscope. According to:

- British Society of Gastroenterology
- HSE Standards and Recommended Practices for Endoscope Reprocessing Units
- European Society of Gastroenterology
- AAMI ST91 American National Standard
- Society of Gastroenterology Nurses and Associates



Wassenburg Medical B.V.  
Edisonring 9  
6669 NA DODEWAARD  
THE NETHERLANDS

T: +31 (0) 488 700 500  
F: +31 (0) 488 453 685  
[info@wassenburgmedical.nl](mailto:info@wassenburgmedical.nl)  
[www.wassenburgmedical.com](http://www.wassenburgmedical.com)