

Specialists in endoscope reprocessing since 1984!



EndoXpert[®] Manager

Data-powered confidence in every cycle

In endoscopy departments, infection prevention and efficiency have the highest priority. Every reprocessing cycle must meet the highest hygiene standards while keeping workflows smooth. Safety starts with infection prevention - powered by smart data management.

Confidence in every cycle

EndoXpert® Manager is an advanced data management platform that registers, organizes and secures all critical process data related to the cycle of endoscope reprocessing within the endoscopy department, storing it in such a way that ensures instant accessibility for analysis and traceability needs. EndoXpert Manager is designed for those who require precision, transparency and compliance. It combines a modern, intuitive interface to streamline workflows and reduce complexity

with advanced security and full compatibility across Wassenburg equipment. With real-time statistics, process graphs, and External Protocol support, EndoXpert Manager provides the data-driven insights required to minimize risks, ensure traceability, and meet international security standards. Built-in proactive maintenance tools and future-ready updates keep the operations efficient and deliver peace of mind.

EndoXpert Manager - when safety, compliance and efficiency matter most

Key features



Compatibility with all Wassenburg® equipment



Modern user interface



Integration with Hospital Information Systems (HIS) and third-party software applications



Digital printer ticket



Continuous improvement with regular updates



Service Tool for configuration



Statistics & Washing - Disinfection process graphs



Enhanced security:

- Input validation
- OAuth support for External Protocol
- Https communication
- Secure login access

EndoXpert® Manager with adaptive licensing

The EndoXpert Manager includes digital printer ticket functionality for seamless documentation. Upgrading to the Advanced license delivers in-depth statistical and graphical insights, enabling deeper analysis of washer-disinfector. For improved process integration, adding

the External Protocol allows EndoXpert Manager to seamlessly connect with hospital information systems or other third-party software applications, supporting efficient workflows and data-driven decision-making.

	EndoXpert Manager	EndoXpert Manager Advanced	Endoxpert Manager Advanced with external protocol
Digital printer ticket	✓	✓	✓
Service Tool	✓	✓	✓
Washer-Disinfector statistics		✓	✓
Washer-Disinfector process graph		✓	✓
Integration with Hospital Information Systems			✓